

Fc-fusion Cell Line Development: Expression and Analytical Strategies

By Leon Pybus¹, Fay Saunders¹, Devika Kalsi¹, Jeffrey Keen¹, Jessica Kane-Fidgeon¹, Daniel Sharkey¹, Yuta Murakami², Tatsuya Matsuura², Kenta Seko²

¹FUJIFILM Diosynth Biotechnologies, Teesside, UK ²FUJIFILM Corporation, Ashigarakami-gun, Kanagawa, Japan

Following the clinical and commercial success of monoclonal antibodies, Fc-based therapies have captured substantial attention in the biopharmaceutical industry; however, their inherent complexity continues to present manufacturability challenges in both upstream and downstream processes. We describe high titer expression of Fc-fusion proteins using FUJIFILM Diosynth Biotechnologies' Apollo[™] X platform coupled with integrated quality assessments, enhancing our support for clients through the predicted growth in Fc-fusion clinical applications.

Key points

Fc-fusion proteins: improving therapeutic potential New developments and future perspectives Manufacturability and quality assessment challenges Apollo[™] X: a versatile expression system for accelerated timelines Platform capabilities Fc-fusion protein productivity assays Clonal expression, purity, and quality assessments

Fc-fusion Proteins: Enhanced Pharmacokinetics and Therapeutic Activity

Fc-fusion proteins are engineered to contain an immunoglobin Fc domain that is directly linked to the protein or peptide of interest (Fig. 1),¹ endowing these hybrid candidates with beneficial pharmacokinetic properties such as increased plasma half-life, stability, avidity, and potency,² and enhancing their therapeutic activity against a wide range of pathologies.^{3,4}





<u>Figure 1</u>. Basic structure of Fc-fusion proteins including the Fc-domain. Adapted from *Duivelshof et al 2021 (J Sep Sci)*.

Although monoclonal antibodies (mAbs) remain the top selling biologics globally,⁵ clinical applications of commercial Fc-based therapies continue to grow (Table 1). In Europe and the US, manufacturing demand for clinical and commercial antibody fusion products is expected to grow at an annual rate of \geq 3% (Table 2), with 30-40% of this demand expected to be outsourced to CDMOs.⁶

<u>Table 1</u>. Fourteen Fc-fusion proteins have been approved in the European Union and United States, and numerous Fc-fusion products are currently in pre-clinical and clinical development.^{1,7}

Product	Company	Fc-fused Protein	2021 sales (\$M)
Eylea	Regeneron Pharmaceuticals Inc.	VEGF Receptors 1 and 2	\$ 9,241.40
Trulicity	Eli Lilly and Co.	Glucagon-like peptide 1 (GLP-1)	\$ 6,174.90
Enbrel	Amgen Inc.	Tumor necrosis factor receptor (TNFR p75)	\$ 5,650.00
Orencia	Bristol Myers Squibb Co.	Cytotoxic T lymphocyte antigen-4 (CTLA-4)	\$ 3,306.00
Eloctate Elocta (EU)	Sanofi SA (Bioverativ Inc.)	Factor VIII	\$ 1,126.80
AlprolIX	Sanofi SA (Bioverativ Inc.)	Factor IX	\$ 694.80
Strensiq	AstraZeneca plc (Alexion)	Catalytic domain of human TNAP	\$ 575.50
Eticovo Benepali (EU)	Samsung Bioepis	Tumor necrosis factor receptor (TNFR p75)	\$ 561.30
Reblozyl	Merck & Co. Inc.	Activin receptor type 2B (ACVR2B)	\$ 551.00
Erelzi	Novartis AG (Sandoz)	Tumor necrosis factor receptor (TNFR p75)	\$ 170.00
Zaltrap	Regeneron Pharmaceuticals Inc.	VEGF Receptors 1 and 2	\$ 91.70
Nulojix	Bristol Myers Squibb Co.	Cytotoxic T lymphocyte antigen-4 (CTLA-4)	\$ 90.00
Arcalyst	Kiniksa Pharmaceuticals Ltd.	IL-1 receptor component, accessory protein	\$ 40.70
Nepexto (EU)	Viatris Inc.	Tumor necrosis factor receptor (TNFR p75)	\$ 50.00



Fc-Fusion Product (1,000 L/Year)	2022	2023	2024	2025	2026
Commercial	261	267	277	290	295
BLA/MAA	0.03	0.1	0.1	0.2	0.3
Phase 3	-	-	-	-	0.02
Phase 2	-	-	-	-	1
Phase 1	-	-	-	-	-
Clinical	13	15	16	17	19
Total	274	282	293	307	315

Table 2. Estimated production growth of Fc-based biologics.

However, despite significant advances in fusion protein technology, progress in bench-to-bedside development of Fc-based therapies remains hindered. Poor physicochemistry can lead to failure of candidates in clinical trials,^{8,9} and optimization of therapeutic outcomes via Fc fragment modification is still an area of intense investigation.⁸ Glycosylation plays a particularly important role in the stability, therapeutic potential, and immunogenicity of mAbs and Fc-fusion products.¹⁰ Assessing glycosylation of Fc-fusion candidates throughout the production process is vital to ensure product stability and proper safety and efficacy profiles.

Fc-fusion manufacturability and quality challenges

The unique structure of Fc-fusion proteins can lead to expression, purity, stability, and quality assessment challenges that must be addressed to meet increasing demand.

Production yields for fusion proteins can be lower than those of mAbs due to less efficient translation or increased susceptibility to proteolytic degradation during production.¹¹ Incomplete or heterogeneous glycosylation can also affect functionality, conformation and/or stability of the fusion product, further complicating manufacturability.⁸ While most therapeutic mAbs contain glycosylation sites in the Fc region and (in some cases) the Fab region, glycosylation of Fc-fusion proteins can also occur in the fusion partners, resulting in diverse and differential glycosylation patterns that have the potential to more significantly impact Fc-fusion product pharmacokinetics, efficacy, and safety profiles.¹¹

The structural diversity and complex glycosylation of Fc-fusion proteins can also present analytical characterization challenges, often requiring product-specific methods in place of available generic or platform technologies that are typically used to assess the quality of biologics.⁴



Apollo™ X: Bridging the Productivity Gap

Expression in mammalian systems under optimal bioreactor^{12,13} and feed supplementation conditions can directly impact Fc-fusion protein conformation, folding, and post-translational modifications, allowing for increased specific productivity, decreased aggregation, and optimal glycosylation.¹⁴

Successful implementation of these strategies requires extensive technical knowledge and practical know-how. FUJIFILM Diosynth Biotechnologies' Apollo™ mammalian Х expression system combines streamlined process development and manufacturability assessment with biomanufacturing-ready cell lines. This allows FUJIFILM Diosynth Biotechnologies to take biopharmaceuticals from pre-clinical to commercial production efficiently and rapidly.

Apollo[™] X was created using an adapted evolution strategy¹⁵ to exploit intrinsic host cell

FEATURES

Host cell line

- DG44-derived CHO host cell line
- Adapted to tailor-made, chemically-defined basal medium and suspension culture
- Superior growth characteristics and expression capability
- Fully characterized cGMP cell bank

Expression Vector

- Double gene vector
- DHFR-based selection system
- Proprietary leader sequence and codon optimization for stable expression and efficient secretion of the target protein

Process

- Robust cloning strategy
- Automated productivity screening and cell isolation remove CLD bottlenecks for high specificity and quality
- Multi-candidate screening and manufacturability assessment package identifies lead candidates and de-risks development process

line heterogeneity and select for superior biomanufacturing potential. A novel workflow aids in the generation of high-quality clonal cell lines with >99% probability of monoclonality in a single round of cloning, balancing speed and quality in cell line development.

The stable and clonal cell lines generated using the Apollo[™] X system are suitable for use in both fed-batch and continuous culture, supporting an intensified DNA-to-research-cell-bank timeline of approximately 10 weeks, and achieving product concentrations of >5 g/L for mAbs (Fig. 2). The process is simple, robust, and scalable for every bioreactor size, and can be used for mAbs and next generation biotherapeutics such as bispecific antibodies and Fc-fusions.



Figure 2. CLD process timeline. Product concentrations of >5 g/L are typically achieved 10 weeks after transfection.

Here, we test the platform's ability to express multiple Fc-fusion products, as well as the quality assessment capabilities of an abbreviated analytical package.

Selecting relevant Fc-fusion proteins for productivity assays

A survey of Fc-fusion manufacturing trends revealed that there are 26 types of Fc-fusion molecules currently approved or in clinical development (Table 3). Analysis of their properties using information contained in the public monoclonal antibodies database IMGT/mAb-DB indicated that 65% of these fusion products contain an IGHG1 Fc fragment. Additionally, these proteins are typically 50-120 kDa in size, are mostly homodimers, and have isoelectric points (pIs) that range between 5-9 (Fig. 3).



Figure 3. Physicochemical properties of Fc-fusion products currently approved or in clinical development.

Based on these findings, we set out to verify the productivity of Fc-fusion molecules that contain an IGHG1 Fc fragment and form bilaterally symmetrical dimer (homodimer) structures in Apollo™ X.

FUJIFUM



<u>Table 3</u>. Characteristics of Fc-fusion-type therapeutics currently approved or in clinical development (Source: Evaluate Ltd 2021).

Number of Products	Fc-fused Protein	Fc Fragment Type	Structure
1	CTLA4	IGHG1	Homo dimer
2	ICOSLG	IGHG1	Homo dimer
15	VEGFR-1 – VEGFR-2	IGHG1	Homo dimer
2	TNFRSF6	IGHG1	Homo dimer
2	TNFRSF13B	IGHG1	Homo dimer
1	CTLA4	IGHG1	Homo dimer
2	peptide 16-mer – peptide 19-mer	IGHG1	Homo dimer
3	F8 – VWF – XTEN	IGHG1	-
1	IL2	IGHG1	Homo dimer*
1	CSF3	IGHG2	Homo dimer
1	-	HyFc	_
1	IL22	IGHG2	Homo dimer
2	IL22	IGHG4	Homo dimer
1	_	—	_
1	GH1 - PEG Linker	IGHG4	Hetero dimer
2	LAG3	IGHG1	Homo dimer
1	coagulation factor IX	IGHG1	Hetero dimer
4	EDA	_	Hexamer
13	TNFRSF1B	IGHG1	Homo dimer
1	SIRPA	IGHG1	Homo dimer
2	_	—	-
1	—	_	-
1	IL1RAP – IL1R1	IGHG1	Homo dimer
2	ACVR2A	IGHG1	Homo dimer
2	peptide 14-mer – peptide 14-mer IGHG1 Homo dim		Homo dimer*
1	TNFRSF1B	IGHG1	Homo dimer
	Number of Products 1 2 15 2 15 2 1 2 1 2 3 1 2 3 1 2 3 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 1 2 2 2 2 2 2	Number of ProductsFc-fused Protein1CTLA42ICOSLG15VEGFR-1 – VEGFR-22TNFRSF62TNFRSF13B1CTLA42peptide 16-mer – peptide 19-mer3F8 – VWF – XTEN1IL21CSF31-1IL222IL221GH1 - PEG Linker2LAG31coagulation factor IX4EDA13SIRPA2.1IL1RAP – IL1R12ACVR2Apeptide 14-mer – peptide 14-mer1TNFRSF1B	Number of ProductsFc-fused ProteinFc Fragment Type1CTLA4IGHG12ICOSLGIGHG115VEGFR-1 – VEGFR-2IGHG12TNFRSF6IGHG12TNFRSF13BIGHG11CTLA4IGHG12peptide 16-mer – peptide 19-merIGHG13F8 – VWF – XTENIGHG11IL2IGHG21IL2IGHG21IL2IGHG41CSF3IGHG21IL2IGHG41Cagalation factor IXIGHG11Coagulation factor IXIGHG11SIRPAIGHG11JIRAP – ILIR1IGHG12ACVR2AIGHG12ACVR2AIGHG11TNFRSF1BIGHG12Peptide 14-mer – peptide 14-merIGHG11TNFRSF1BIGHG1

*N-terminus Fc



To assess Apollo[™] X expression capabilities and evaluate the effects of MW (molecular weight) and pl (isoelectric point) on Fc-fusion protein productivity, we selected seven molecules with varying MW and pl based on amino acid sequence information (Fig. 4, Table 4 below). Productivity results should therefore be representative of a wide range of Fc-fusion proteins that are currently approved for clinical use or still in development.



Figure 4. Fc-fusions selected for analysis.

Assessing Fc-fusion productivity

Stable expression cell lines were established using a

single gene vector (SGV) for each Fc-fusion product. SGVs were linearized and transfected into Apollo[™] X cells via electroporation at uniform amounts to minimize the effect of gene insertion number on productivity. Stable expression pools were generated by selection with MTX and expanded into 125 mL volume non-baffled vented flasks for fed-batch culture. Standard methods were used to assess cell viability, Fc-fusion product concentration, and purification results.

Cell culture was performed using Fc-fusion transfectant pools to assess productivity in the ApolloTM X platform. Six of the seven Fc-fusion products tested reached productivity levels of 1 g/L or more using fed-batch culture of transfectant pools (Table 4); titer is expected to be higher in clonal cell lines. The viable cell density (VCD) peak of each Fc-fusion expressing cell pool was $^{40} \times 10^{6}$ cells/mL, similar to those attained in mAb production.

Fc-fusion	Theoretical MW [kDa]	рІ	IVCD [x10 ⁶ (day · cells)/mL]	Titer* [g/L]
Fc01	79	5.7	319	1.4
Fc02	97	8.2	312	1.1
Fc03	71	8.1	306	1.2
Fc04	64	7.2	358	2.2
Fc05	102	7.2	274	1.6
Fc06	78	6.0	293	1.6
Fc07	201	6.9	330	0.3

Table 4. Fed-batch culture of Fc-fusion proteins selected for productivity assays.

*Titer was measured with CedexBio and corrected by theoretical MW



Examining the effect of Fc-fusion properties on productivity

A preliminary study assessing expression of Fc-fusion proteins of varying size suggested that MW has a significant impact on Fc-fusion productivity. On the other hand, pl was hypothesized to influence degradation and aggregation during cell culture and downstream purification processes. We assessed the effects of MW, pl, and mRNA levels on production of the Fc-fusion proteins selected (Fig. 5).



<u>Figure 5</u>. mRNA levels, MW, pl, and Fc-fusion product titer in transfectant pools.

MW and mRNA levels were closely correlated with titer. These results were not unexpected: transcription has considered long been а dominant factor in protein expression control, and correlations between transcript levels and mAb productivity have been observed.16,17,18,19

In general, Fc-fusion molecules of \leq 110 kDa achieved productivities of 1 g/L or more in ApolloTM X transfectant pools. On the other hand, we observed no effect of pl on

productivity (Fig. 5D), suggesting that pl-mediated aggregation is not an issue for these Fc-fusions when the Apollo™ X workflow is used, or that it is not at play within the range tested.

Fc02: Clonal cell line development and product analysis

To further explore the Fc-fusion expression capabilities of Apollo[™] X and the analytical power of our abbreviated package, we took a closer look at Fc02 productivity and quality assessment. Fc02's inherent structural complexity and posttranslational modifications make it a good candidate to map the possible challenges we may encounter in the production and analysis of Fc-fusion proteins.



Cell line development

48 cell lines expressing Fc02 were generated using the Apollo[™] X platform's cell line development (CLD) process and evaluated in the Ambr®15 fed-batch screen over a 14-day period.

Final product concentrations ranged from 3.4 to 4.3 g/L for the 10 highest producers (Fig. 6), suggesting that under clonal conditions, the platform's standard development workflow is sufficient to produce Fc-fusion proteins at titers comparable to those of mAbs.



<u>Figure 6</u>. Fc02 CLD results. Protein concentration was determined using Protein A UPLC.

Protein expression and purity assays

MW and pI are often used to estimate the purity of antibody preparations²⁰; however, the increased variety and complexity of Fc fusion glycosylations can greatly affect theoretical pI and MW calculations.^{21,22} CE-SDS analysis of Fc02 before and after deglycosylation by enzymatic digestion revealed that native Fc02 exhibits significant deviation from the theoretical MW, whereas deglycosylation resulted in a more accurate approximation with a narrower peak profile (Fig. 7).



The increased variation observed in the intact form is suggestive of heavy glycosylation, so we examined glycan subtypes in more detail.

Figure 7.

Approximation of theoretical MW by CE-SDS varies as a function of glycosylation. [Fc02 theoretical MW = 97 kDa (dimer)].



Assessing Fc-fusion product by mass spectrometry

Fc02 product was evaluated for intact mass as detailed in Table 5a. Further investigations into glycosylation were performed using the methodology in Table 5b.

Intact mass analysis of Fc02 supports sequence integrity with observed masses agreeing well with expected mass (Table 6) while antennary type 2 glycans were the most predominant forms, as typically observed with mAbs (Table 7). Mannose, sialic acid, fucose, and galactose modifications (Fig. 8) have significant effects on IgG-based therapeutics SO monitoring glycans is crucial.²² Product quality assays such as mass spectrometry are needed during development process and manufacturing avoid to irreproducibility misleading or observations in upstream and downstream applications. Here, we



Figure 8. Common N-linked glycan types in mAb and Fc-fusion glycosylation.

show that FUJIFILM Diosynth Biotechnologies has the capability to conduct this type of analysis in flexible Fc-fusion process development programs.



Table 5a. Intact mass methodology.

BioAccord System:	Waters Acquity UPLC I-Class with TUV detector coupled to Acquity RDa Mass Spectrometer; Spectral data processed using Waters UNIFI software.		
Analytical column:	Acquity UPLC Protein BEH C4, 300 Å, 1.7 μm, 2.1x50 mm		
Mobile phase:	MPA: 0.1% Formic Acid in LC-MS UHQH ₂ O MPB: 0.1% Formic Acid in LC-MS Acetonitrile		
Sample preparation:	Intact: Dilute to 1 mg/mL in UHQH ₂ O Deglycosylation: Diluted to 1 mg/mL treated overnight with PNGase F Reduced: Intact/deglycosylated sample treated with guanidine hydrochloride to denature and reducing agent		
Sample loading:	2 μL		
Flow rate:	0.2 mL/min		
Run time:	10 minutes		
Autosampler temperature:	6°C		
Column temperature:	80°C		
UV detection:	215 nm		

Table 5b. Glycan analysis methodology.

BioAccord System:	Waters Acquity UPLC I-Class with fluorescence detector coupled to Acquity RDa Mass Spectrometer; Spectral data processed using Waters UNIFI software.		
Analytical column:	Acquity Premier Glycan BEH Amide, 150 mm x 2.1 mm, 1.7 μ m		
Mobile phase:	Mobile Phase A: 50 mM Ammonium Formate pH 4.4 Mobile Phase B: Acetonitrile		
Sample preparation:	Samples were prepared using the Waters RapiFluorMS standard protocol		
Sample loading:	25 μL		
Flow rate:	0.4 mL/min and 0.2 mL/min		
Run time:	55 minutes		
Autosampler temperature:	10°C		
Column temperature:	60°C		
Fluorescence detection:	Excitation=265 nm, Emission=425 nm		



Table 6. Intact mass analysis of deglycosylated Fc02.

Protein name	Expected mass (Da)	Observed mass (Da)	Mass difference (Da)
Fc02	96897.2	96909.1	+11.9
Fc02 (reduced)	48458.7	48461.4	+2.7

Table7. Fc02 glycosylation analysis.

Species	% Total
Total Sialylated Glycans	11%
Mono Sialylated	9%
Di Sialylated	2%
Non-Sialylated Glycans	89%
Defucosylated Glycans	39%
Antennary 1 Glycans	10%
Antennary 2 Glycans	56%
Antennary 3 Glycans	3%
High Mannose Glycans	29%

Perspectives

Fc-fusion proteins have emerged as next generation biotherapeutics due to their improved pharmacokinetic properties and therapeutic activity, but despite intense investigation, development remains a challenge. The inherent complexity and structural diversity of these molecules can result in low expression, molecule instability, and inefficient processes. In this study, we evaluated the adaptability of the Apollo[™] X mammalian expression system for expression of seven IGHG1-type Fc-fusion proteins of varying physicochemical characteristics.

Fed-batch culture of transfectant pools resulted in productivity levels ≥ 1 g/L for Fc-fusion proteins with a theoretical MW ≤ 110 kDa, a characteristic which describes most of the Fc-based therapeutics currently approved or in development. We believe that productivity of Fc-fusions with higher molecular weight may be improved by selecting clones with high Gene of Interest (GOI) mRNA expression or by taking measures to stabilize mRNA levels.



These data suggest that (i) mRNA transcript levels could be used to screen for cell lines with increased Fc-fusion productivity, and (ii) vector²³, cell line²⁴ and/or process optimization²⁵ strategies to increase mRNA abundance via increased transcription or mRNA stability. These strategies are options for future developments within FUJIFILM Diosynth Biotechnologies.

Having the ability to evaluate Fc-fusion product quality quickly and reliably is crucial to ensure safety and efficacy in downstream applications. Antibody preparation purity is often assessed via MW and pI estimations; however, these conventional methods can be rendered impractical given the extensive glycosylation of Fc-fusion proteins. Glycosylation has implications on Fc-fusion stability, pharmacokinetics, and immunogenicity, and must be closely monitored throughout biologic manufacturing.²⁶ Here, we confirm the effect of glycosylation on MW estimations, with CE-SDS indicating a significant shift in mass for deglycosylated product in comparison to intact product and mass spectrometry of deglycosylated product confirming expected mass. We propose an abbreviated workflow that allows for accurate and reliable evaluation of product homogeneity and quality.

The versatile Apollo[™] X mammalian expression system combines cell line development, automated productivity screening, and manufacturability assessment services into one effective, integrated solution for monoclonal and bispecific antibodies, and Fc-fusion proteins. Our ability to adapt well established methods to assess product quality early in the process chain further streamlines Fc-fusion protein development, propelling these products through critical milestones with quality and patient safety at the fore.

Authors

Dr. Kalsi is a Principal Scientist in Mammalian Upstream Process Development, Dr. Saunders is the Head of Mammalian Cell Culture, Dr. Pybus is a Senior Staff Scientist in Mammalian Upstream Process Development, Jessica Kane-Fidgeon is a Principal Scientist in Analytical Development, Daniel Sharkey is a Research Scientist in Analytical Development, Dr. Keen is a Senior Staff Scientist in Analytical Development at FUJIFILM Diosynth Biotechnologies, Belasis Avenue, Billingham, TS23 1LH, UK.

Mr. Murakami is a Research Manager in Cell Science and Engineering,

Mr. Matsuura is a Senior Scientist in Mammalian Cell Line Development,

Mr. Seko is a Scientist in Mammalian Cell Line Development

at FUJIFILM Bio Science and Engineering Laboratories, Ashigarakami-gun, Kanagawa, Japan.

References

¹ Duivelshof BL, Murisier A, Camperi J, Fekete S, Beck A, Guillarme D, D'Atri V. 2021. Therapeutic Fc-fusion proteins: Current analytical strategies. *J Sep Sci* 44(1):35-62, PMID: 32914936,

https://doi.org/10.1002/jssc.202000765.

² Czajkowsky DM, Hu J, Shao Z, Pleass RJ. 2012. Fc-fusion proteins: new developments and future perspectives. *EMBO Mol Med* 4(10):1015-28, PMID: 22837174, https://doi.org/10.1002/emmm.201201379.

³ Roopenian DC, Akilesh S. 2007. FcRn: the neonatal Fc receptor comes of age. *Nat Rev Immunol* 7(9):715-25, PMID: 17703228, https://doi.org/10.1038/nri2155.

⁴ Kontermann RE. 2011. Strategies for extended serum half-life of protein therapeutics. *Curr Opin Biotechnol* 22(6):868-76, PMID: 21862310, https://doi.org/10.1016/j.copbio.2011.06.012.

⁵ Lu RM, Hwang YC, Liu IJ, Lee CC, Tsai HZ, Li HJ, Wu HC. 2020. Development of therapeutic antibodies for the treatment of diseases. *J Biomed Sci* 27(1):1, PMID: 31894001, https://doi.org/10.1186/s12929-019-0592z.

⁶ bioTRAK Database. 2021. SUPPLY AND DEMAND TRENDS: MAMMALIAN BIOMANUFACTURING INDUSTRY OVERVIEW. [Website]. Boston, MA, USA: BioProcess Technology Group (BDO).

https://www.bdo.com/getattachment/5091b97d-f052-4951-ab84-d88e46244b0a/attachment.aspx?BPI-Mammalian-Biomanufacturing-Industry-Overview.pdf. [Accessed 20 Aug 2022].

⁷ Jafari R, Zolbanin NM, Rafatpanah H, Majidi J, Kazemi T. 2017. Fc-fusion proteins in therapy: An updated view. *Curr Med Chem* 24(12):1228-1237, PMID: 28088904, https://doi.org/10.2174/0929867324666170113112 759.

⁸ Yang C, Gao X, Gong R. 2018. Engineering of Fc fragments with optimized physicochemical properties implying improvement of clinical potentials for Fc-based therapeutics. *Front Immunol* 8:1860, PMID: 29375551, https://doi.org/10.3389/fimmu.2017.01860.

⁹ Rath T, Baker K, Dumont JA, Peters RT, Jiang H, et al. 2015. Fc-fusion proteins and FcRn: structural insights for longer-lasting and more effective therapeutics. *Crit Rev Biotechnol* 35(2):235-54, PMID: 24156398, https://doi.org/10.3109/07388551.2013.834293.

¹⁰ Liu L. Antibody glycosylation and its impact on the pharmacokinetics and pharmacodynamics of monoclonal antibodies and Fc-fusion proteins. *J Pharm Sci* 104(6):1866-1884, PMID: 25872915, https://doi.org/10.1002/jps.24444.

¹¹ Challener CA. 2017. Fusion proteins pose manufacturability challenges. *BioPharm International* 30(5):30,

https://www.biopharminternational.com/view/fusionproteins-pose-manufacturability-challenges.

¹² Trummer E, Fauland K, Seidinger S, Schriebl K, Lattenmayer C, et al. 2006. Process parameter shifting: Part II. Biphasic cultivation-A tool for enhancing the volumetric productivity of batch processes using Epo-Fc expressing CHO cells. *Biotechnol bioengin* 94(6): 1045-52, PMID: 16736532, https://doi.org/10.1002/bit.20958.

¹³ Estes B, Hsu YR, Tam LT, Sheng J, Stevens J, Haldankar R. 2015. Uncovering methods for the prevention of protein aggregation and improvement of product quality in a transient expression system. *Biotechnol Prog* 31(1): 258-67, PMID: 25395220, https://doi.org10.1002/btpr.2021.

¹⁴ Jazayeri JA, Carroll GJ. 2008. Fc-based cytokines: prospects for engineering superior therapeutics. *BioDrugs* 22(1):11-26, PMID: 18215087, https://doi.org/10.2165/00063030-200822010-00002.

¹⁵ Majors BS, Chiang GG, Betenbaugh MJ. 2009. Protein and genome evolution in Mammalian cells for biotechnology applications. *Mol Biotechnol* 42(2):216-23, PMID: 19367473, https://doi.org/10.1007/s12033-009-9156-x.

¹⁶ Gutierrez JM, Feizi A, Li S, Kallehauge TB, Hefzi H, et al. 2020. Genome-scale reconstructions of the mammalian secretory pathway predict metabolic costs and limitations of protein secretion. *Nat Commun* 11(1):68, PMID: 31896772, https://doi.org/10.1038/s41467-019-13867-y.

¹⁷ Jiang Z, Huang Y, Sharfstein ST. 2006. Regulation of recombinant monoclonal antibody production in chinese hamster ovary cells: a comparative study of gene copy number, mRNA level, and protein expression. *Biotechnol Prog* 22(1):313-8, PMID: 16454525, https://doi.org/10.1021/bp0501524.

¹⁸ Dorai H, Csirke B, Scallon B, Ganguly S. 2006. Correlation of heavy and light chain mRNA copy numbers to antibody productivity in mouse myeloma production cell lines. *Hybridoma* (Larchmt) 25(1):1-9, PMID: 16475875, https://doi.org/10.1089/hyb.2006.25.1.

¹⁹ Lee CJ, Seth G, Tsukuda J, Hamilton RW. 2009. A clone screening method using mRNA levels to determine specific productivity and product quality for monoclonal antibodies. *Biotechnol Bioeng* 102(4):1107-18, PMID: 18985612, https://doi.org/10.1002/bit.22126.

²⁰ Raynal B, Lenormand P, Baron B, Hoos S, England P. 2014. Quality assessment and optimization of purified protein samples: why and how? *Microb Cell Fact* 13:180, PMCID: PMC4299812, https://doi.org/10.1186/s12934-014-0180-6.

²¹ Wang G, de Jong RN, van den Bremer ETJ, Parren PWHI, Heck AJR. 2017. Enhancing Accuracy in Molecular Weight Determination of Highly Heterogeneously Glycosylated Proteins by Native Tandem Mass Spectrometry. *Anal Chem* 89(9):4793-4797, PMID: 28383250, https://doi.org/10.1021/acs.analchem.6b05129.

²² Boune S, Hu P, Epstein AL, Khawli LA. 2020. Principles of *N*-Linked Glycosylation Variations of IgG-Based Therapeutics: Pharmacokinetic and Functional Considerations. *Antibodies* (Basel) 9(2):22, PMID: 32532067, https://doi.org/10.3390/antib9020022.



 ²³ Hung F, Deng L, Ravnikar P, Condon R, Li B, et al. 2010. mRNA stability and antibody production in CHO cells: improvement through gene optimization. *Biotechnol J* 5(4):393-401, PMID: 20222103, https://doi.org/10.1002/biot.200900192.

²⁴ Brown AJ, James DC. 2016. Precision control of recombinant gene transcription for CHO cell synthetic biology. *Biotechnol Adv* 4(5):492-503, PMID: 26721629, https://doi.org/10.1016/j.biotechadv.2015.12.012. ²⁵ Gutiérrez-González M, Latorre Y, Zúñiga R, Aguillón JC, Molina MC, Altamirano C. 2019. Transcription factor engineering in CHO cells for recombinant protein production. *Crit Rev Biotechnol* 39(5):665-679, PMID: 31030575,

https://doi.org/10.1080/07388551.2019.1605496.

²⁶ Edwards E, Livanos M, Krueger A, Dell A, Haslam SM, et al. 2022. Strategies to control therapeutic antibody glycosylation during bioprocessing: Synthesis and separation. *Biotechnol Bioeng* 119(6):1343-1358, PMID: 35182428, https://doi.org/10.1002/bit.28066.